

# ETSimulations User Manual

Kyung Min Shin

Ver 0.2.2 - March 2020

## Contents

<b>1 Getting started</b>	<b>1</b>
External dependencies . . . . .	1
Setting up ETSimulations . . . . .	1
Running a simulation . . . . .	1
Processing a simulated dataset . . . . .	2
<b>2 Simulation configuration parameters</b>	<b>2</b>
<b>3 Dataset generation outputs</b>	<b>3</b>
<b>4 The EMAN2 Processor</b>	<b>4</b>
steps_to_run . . . . .	5
particles_coordinates_file . . . . .	5
unbinned_boxsize . . . . .	5
Parameters objects . . . . .	5
Running the generated script . . . . .	6
<b>5 The IMOD Processor</b>	<b>6</b>

# 1 Getting started

The ETSimulations package is designed to be a tool to facilitate generation and processing of simulated cryo-electron tomography data. As of this preliminary version, only the generation of simulated tilt stacks is well-supported and covered by this manual.

## External dependencies

The following external software is used by ETSimulations and should thus be installed first:

- Chimera - <https://www.cgl.ucsf.edu/chimera/>
- TEM-simulator - <http://tem-simulator.sourceforge.net/>
- Python 3.3 or greater

## Setting up ETSimulations

To begin, navigate to a directory in which you want to install the package and run the following commands:

```
git clone https://github.com/kmshin1397/ETSimulations.git
```

```
cd ETSimulations
```

```
python3 -m venv env
```

```
source env/bin/activate
```

```
pip install -r requirements.txt
```

This will have downloaded the ETSimulations code and set up the necessary Python packages and environment.

## Running a simulation

Most all general parameter set up is done through a YAML file which is passed in as an argument to the main **ets\_generate\_data.py** program. An example such YAML file is provided in the ETSimulations directory as **configs.yaml**.

Detailed set-up with regards to the characteristics of the particles simulated are controlled through custom "**Assembler**" Python classes which can define a series of Chimera commands to open, manipulate, and combine one or more source maps, i.e. various proteins, into a fake particle source. (This is done for the T4SS simulations in the `src/assemblers/t4ss_assembler.py` file) Custom arguments to pass

along to your custom particle Assemblers can also be defined in the configuration YAML file, as shown in the example.

To actually run a set of simulations, run (assuming you've activated the virtual environment as indicated above - the **source ...** command):

```
python src/ets_generate_data.py -i configs.yaml
```

More details on the parameters available can be found in the next section.

## Processing a simulated dataset

Another program called **ets\_process\_data.py** is provided to facilitate processing of the raw datasets generated by the simulations. Here as well, we provide input via a configuration YAML file. This time parameters are given to modules called "**Processors**". Each Processor module implements a function to take the inputted configurations and automatically set up tomogram reconstruction and sub-tomogram averaging. For example, the EMAN2 processor (the only one implemented so far) will automatically create an EMAN2 project and generate a Python script which will run all desired EMAN2 processing steps (i.e. `e2import.py`, `e2tomogram.py`, etc.) in order.

The YAML file provided to the `ets_process_data.py` run should have:

- **root**: The path to the folder which contains the `raw_data` and `sim_metadata.json` files created by `ets_generate_data.py`
- **name**: The project name - should be the same one used for the simulations
- **processors**: A list of processor objects which have a "name" and "args" arguments. The specifics to what goes into these fields can be found in more details in the section of the manual particular to that processor.

To run the processor set up, run:

```
python src/ets_process_data.py -i processor_configs.yaml
```

## 2 Simulation configuration parameters

- **tem\_simulator\_executable** : The file path to the TEM-simulator executable
- **chimera\_exec\_path** : The file path to your Chimera installation
- **model** : The path to the main particle source file (Doesn't actually matter for T4SS simulations because I'm bypassing this and putting together a bunch of different source maps)

- **root:** The project root directory in which to generate simulations
- **config:** The TEM-simulator configuration text file to apply to each simulation. An example is provided in the templates folder.
- **coord:** The TEM-simulator particle coordinates text file to use as a reference for placing particles in each generated tiltseries. An example for this is also provided in the templates folder.
- **num\_stacks:** The number of tilt stacks to generate
- **name:** A name for the project
- **num\_cores:** The number of parallel cores to utilize
- **apix:** The pixel size to give to generated stacks, in nm
- **num\_chimera\_windows:** The number of Chimera instances to spawn to drive particle assembly. Creating more windows will clutter your display more, but can make simulations run faster. If your particle assembly does not use a lot of Chimera commands/spend a lot of time running Chimera commands, then having multiple Chimera windows may not be necessary
- **defocus\_values:** A list of defocus values to assign to simulations. The values will be evenly assigned across the dataset, i.e. if 2 defocus values are provided half of the simulated stacks in a dataset will be given one and the other given the other.
- **bead\_map :** The MRC map representing fake gold beads to scatter throughout the tilt stacks. An example is provided in the templates folder
- **email :** An email address to send completion notifications to

### 3 Dataset generation outputs

Running the **ets\_generate\_data.py** program as discussed in section 1.3 will result in a **raw\_data** folder being created in the project directory specified in the configurations. In the **raw\_data** folder, each tiltseries will get its own sub-directory titled {name}\_{stack number}. In each sub-directory, you will find a no-noise version of the stack and a normal noisy version.

The other important output to note is the `sim_metadata.json` file. This is a JSON file containing metadata for each tiltseries generated, including custom metadata that can be saved from your custom Assembler. For example, the T4SS Assembler saves the random orientations and random shifts/angles away from the centered/perpendicular positions for each component of the simulated particle which were generated during the run. An easy way to interact with and retrieve this information is the Python json module which can load this json as a Python dictionary, i.e.:

```
import json
metadata = json.load(open("sim_metadata.json", "r"))
```

## 4 The EMAN2 Processor

The EMAN2 Processor, found in `processors/eman2_processor.py`, is implemented to automatically set up cryo-ET processing of generated simulated data. It works by taking a pre-written Python template script, which lists the EMAN2 programs needed to run to take a series of raw tilt stacks to a subtomogram average, and filling in the desired parameters to those programs based on configurations provided in the YAML input file. This template script can be found at `templates/eman2/eman2_process.py`. Any argument which can be given to the EMAN2 programs called can be given in the YAML file as part of the parameters object for the program. Note the example configurations below.

```
processors: [  
  {name: "eman2",  
    args:  
    {  
      particle_coordinates_file: "/Users/kshin/Documents/repositories/ETSimulations  
/templates/eman2/T4SS_coords_3by3.txt",  
      unbinned_boxsize: 128,  
      steps_to_run: ["import", "reconstruct"],  
      e2import_parameters : {  
        "import_tiltseries": "enable",  
        "importation": "copy",  
        "apix": 2.83,  
        "boxsize": 64  
      },  
      e2tomogram_parameters : {  
        "tltstep": 3,  
        "tltax": -90,  
        "npk": 10,  
        "tltkeep": 0.9,  
        "outsize": "1k",  
        "niter": "2,1,1,1",  
        "pkkeep": 0.9,  
        "bxsx": 64,  
        "pk_mindist": 0.125,  
        "filterto": 0.45,  
        "rmbeadthr": 10.0,  
        "threads": 48,  
        "clipz": 350,  
        "notmp": "enable"  
      }  
    }  
  ]  
}]
```

Like the example above, you can include EMAN2 processing in your processing run by adding a pro-

cessor argument to the "processors" list with the name "eman2" and suitable args.

## **steps\_to\_run**

In the "args" field, the "steps\_to\_run" lists the processing steps which will be run when you eventually run the script generated by the processor (Note that the functions will still exist in the script to run steps not included in the list, they just won't be executed unless you manually enable them back). The available steps taken by the steps\_to\_run are "import", "reconstruct", "estimate\_ctf", "extract", "build\_set", "generate\_initial\_model", and "3d\_refinement". Note also that the order matters in this list (i.e. reconstructing before importing tiltseries will result errors).

## **particles\_coordinates\_file**

The particles\_coordinates\_file argument is used to define a text file containing the 3D coordinates of the particles within the simulated tomograms. This will be used during the "extract" processing step to transfer particle coordinates to the EMAN2 project as if picked in Boxer so that the sub-volumes can be extracted for sub-tomogram averaging. The coordinates in this file should be in pixels, based on the EMAN2 conventions (with the origin being in the center of the axis instead of the lower left, for example). An example such file is provided at templates/eman2/T4SS\_coords\_3by3.txt.

## **unbinned\_boxsize**

This defines the boxsize to record for particles when importing the particles from the particles coordinates file into the EMAN2 tomogram info files. This should be in terms of unbinned pixels, i.e. the scale matching the original tiltseries imported.

## **Parameters objects**

Each EMAN2 program is given its own \*\_parameters section in the "args" field, listing all the command line arguments that would be passed in if calling these programs normally. For example, arguments to e2import.py would be listed in the e2import\_parameters field as shown. Arguments which are just flags instead of taking a value, such as the "-help" option available in all these programs, should be put in to the configuration section with the special value of "enable" as can be seen in the example above. The full list of EMAN2 programs exposed by the EMAN2 Processor, and thus able to take their own \*\_parameters section is:

- e2import.py : e2import\_parameters
- e2tomogram.py : e2tomogram\_parameters
- e2spt\_tomoctf.py : e2spt\_tomoctf\_parameters
- e2spt\_extract.py : e2spt\_extract\_parameters
- e2spt\_buildsets.py : e2spt\_buildsets\_parameters
- e2spt\_sgd.py : e2spt\_sgd\_parameters
- e2spt\_refine.py : e2spt\_refine\_parameters

## Running the generated script

The generated EMAN2 processing script that is outputted by `ets_process_data.py` will be located in the newly created EMAN2 project directory in the `processed_data` folder created. This will be a normal Python script you can run, albeit requiring Python 3, using: **`python3 eman2_process.py`**. Note that Python 3 is only used for the proper file IO and kicking off EMAN2 programs. The EMAN2 programs themselves will be run using Python 2 as EMAN2 is still using Python 2 officially.

Another important thing to note is that the created `eman2_process.py` script is not meant to be a rigid program. It has been designed to be easily modifiable - all parameters originally passed in are located towards the top of the file. Thus, the script can be easily opened and edited as necessary, such as if the `steps_to_process` should be modified to pick up from where an error interrupted the processing. For additional clarity and potential modification, a simple text file containing the raw command-line versions of the EMAN2 commands handled by the `eman2_process.py` will be created as well in a file called `eman2_process_commands.txt`.

## 5 The IMOD Processor

The IMOD Processor, found in `processors/imod_processor.py`, is implemented to facilitate processing of the generated simulation dataset with the IMOD software. Specifically, a project directory is created and set up, along with necessary files, for running the **batchruntomo** program (<https://bio3d.colorado.edu/imod/doc/man/batchruntomo.html>) for processing the simulated stacks in batch. The IMOD Processor as well will take its inputs from the configuration YAML passed into `ets_process_data.py`. Consider the following example:

```
processors: [
  {
    name: "imod",
    args: {
      start_step: 0.0,
      end_step: 8.0,
      num_fiducials: 10,
      tilt_axis: -90,
      apix: 0.283,
      fiducial_method: "autofidseed"
    }
  }
]
```

The IMOD Processor, like all others, have the **name** argument ("imod") and an **args** object filled with parameters. The parameters exposed here are abstracted versions of some of the more common ones when processing with batchruntomo. Specifically, we have:

- **start\_step**: The batchruntomo processing start step (the list of steps can be found on the IMOD batchruntomo documentation page)

- `end_step`: The batchruntomomo processing end step (the list of steps can be found on the IMOD batchruntomomo documentation page)
- `num_fiducials`: The number of fiducials to track with either Raptor or the automatic fiducial seeding/tracking for alignment
- `tilt_axis`: The tiltseries tilt axis you specified in the TEM-Simulator configs
- `apix`: The pixel size (in nm)
- `fiducial_method`: Either "autofidseed" or "raptor", determining how beads are seeded and tracked for alignment
- `custom_template`: (Note this is not included in the example above as putting this option in overwrites the others) A filepath to a batchruntomomo template file (.adoc) to be used instead of the default one generated with the other options. This can be used to control any parameters open to batchruntomomo, but will probably require more experimentation/advanced knowledge to get it to work.

An important thing to note is that stacks generated by the simulations will most likely not have enough signal to have the IMOD coarse alignment with cross-correlation step work. Instead, the cross-correlation is likelier to shift tilts extremely out of alignment than to make any improvements. Thus, the IMOD Processor is designed to automatically skip the coarse alignment (steps 2.0 and 3.0) and fake having done it by moving in fake versions of its outputs like the .rawtilt file. It is possible to override this, by first doing a run with starting and ending steps 0.0 to 1.0, and then doing a second run with starting step 2.0.

Additionally, as mentioned, fake versions of the rough alignment are used in place of doing the cross-correlation. These fake files are located in the templates/imod folder and can be edited if desired. For example, the rawtilt file is for a 2-degree increment, -54 to 54 degrees tilt scheme and should be edited if using a different tilt scheme with your TEM-Simulator.